


A different approach to hemostasis.



NOW available in the USA.

Hemospray[®]
ENDOSCOPIC HEMOSTAT

For more information on Hemospray, visit:
hemospray.cookmedical.com



Images courtesy of Dr. John Morris,
Glasgow Royal Infirmary, Glasgow Scotland.



Hemospray®

ENDOSCOPIC HEMOSTAT

This device is used for hemostasis of nonvariceal gastrointestinal bleeding. Supplied sterile and is disposable – intended single use only.

- **Nonthermal:** No tissue changes occur as sometimes experienced with thermal modalities.
- **Nontraumatic:** Since no force is applied at the treatment site, the powder minimizes the risk of tissue trauma as experienced with other modalities.
- **Noncontact:** Aerosol delivery system eliminates the need for direct mechanical or contact forces with the treatment site by the delivery system.
- **Nonspecific targeting:** Powder is sprayed toward the source of the bleed, so it does not require the en face approach and the precise targeting of other modalities.

Order Number	Reference Part Number	Catheter Length cm	Catheter Size Fr	Minimum Accessory Channel mm
G56572	HEMO-7	220	7	2.8
G21049	HEMO-10	220	10	3.7

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Summary of Clinical Data: Clinical data summary information that was, in part, the basis for granting the de novo can be found on the Cook Medical website at [CookMedical.com/HemosprayData](https://www.cookmedical.com/HemosprayData)

Customer Service

EU Website: [cookmedical.eu](https://www.cookmedical.eu)
 EDI: [cookmedical.eu/edi](https://www.cookmedical.eu/edi)
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: [cookmedical.com](https://www.cookmedical.com)
 EDI: [cookmedical.com/edi](https://www.cookmedical.com/edi)
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4